

SuperSTEP™

Nitinol Bone Fixation System

STABLE

High Strength and force insure
stable fixation

SIMPLE

ZeroPrep™ pre-loaded insertion system
with instantly active implant

VERSATILE

Sizes for mid and rear foot
reconstructions

FUSION


Dynamic compression for
bone healing



METRIC®
BRINGING BONES TOGETHER™

SuperSTEP™

Design Features



Patented ZetaArc™

- Unique variable step height
- High strength and compression forces
- Minimally invasive push in technique
- Shortening bridge and inward moving leg
- Most uniform compression available
- Low bone profile
- Various STEP heights: 6mm, 8mm, 10mm

P/N 20063 STEP-06

Dimensions: STEP-06 20 x 20 x 20 2x3mm

Compression Force: 21.2 lbf

Reduction: Bridge 17% Legs 42%

P/N 20083 STEP-08

Dimensions: STEP-08 20 x 20 x 22 2x3mm

Compression Force: 21.2 lbf

Reduction: Bridge 17% Legs 42%

P/N 20103 STEP-10

Dimensions: STEP-10 20 x 20 x 22 2x3mm

Compression Force: 21.2 lbf

Reduction: Bridge 17% Legs 42%

ZeroPrep™ Pre-loaded Insertion System

- Included disposable instrumentation
- Room temperature storage
- Instantaneously active
- Enhances infection control
- Worker safety considerations
- Logistically simple and ready



**We bring you superior fixation,
Learn more at MetricMD.com**

© 2020 Metric Medical Devices, Inc. 846 Silver Springs, Helotes, TX 78023
(830) 535-6300 A03-008-01 IR

U.S. Federal Law restricts this device to sale by or on the order of a physician
Patents 8,721,646 B2; 9, 451,955 B2; 9,743,926; 9,451,957; 9,339,268; 9,017,331;
10,512,549; E.P. 2,736,41 B1; 10,448,979, other patents pending